



Complete Medical Australasia

ORDER DETAILS

Order Date

Customer Name

Phone

Email

Address

Suburb/City

Email

ORDER REFERENCE NO. (OFFICE USE)

Date Required

Fax

State

Postcode

SCRIPT DETAILS

Don't have a script? Simply call your eyecare professional for your script ensuring you ask for your Pupillary Distance (PD).
Alternatively, send us a spare pair of your own glasses and we will work it out for you.

	Sphere	Cylinder	Axis	Prism	Add	Pupillary Distance Distant	Near	Fitting Heights at PD above rim
R								R
L								L

LENS TYPE ☐ Single Vision ☐ Bifocal ☐ Lenses without Power (Plano) ☐ Progressive

FRAME DETAILS ☐ Own Frame ☐ Rx Safety Stock Frame

ADDITIONAL COMMENTS / INFORMATION

IMPORTANT INFORMATION

- Leaded glass has a SG of 5.1g/cc and lenses are heavier than normal prescription lenses
- Frames are available from our website, or customers can supply their own frame (subject to inspection)
- Own frames must not be nylon or rimless
- We reserve the right to inspect the frame prior to accepting job as some frames such as false front frames are unsuitable and can be difficult to identify
- Wrap frames are unsuitable for scripting with leaded glass, other than unpowered plano lenses
- Leaded glass lenses are unsuitable for multi-coating and chemical hardening processes

SCRIPT RANGES

	Sph Power	Cyl Power	Add Power	
Single Vision	Plus	+6.00D	-4.50D max	N/A
	Minus	-9.00 D net power including cylinder power		
Bifocal D28	Plus	+6.00D	-4.50D max	+1.00D to +3.00D
	Minus	-9.00 D net power including cylinder power		+1.00D to +3.00D
Progressive	Plus	-4.00D to +5.00D	-4.50D max	+1.00D to +3.00D
	Minus	* To be Advised		+1.00D to +3.00D

- Progressive blanks are made to order and take 14 to 21 days
- Progressive and D28 lenses are laminated and are thicker as a consequence
- Powers outside script range specified only after negotiation and depends upon frame size
- Minimum lens thickness: 3mm lead glass

CUSTOMER SIGNATURE

Please complete all sections and send this form along with your original script to:

PH: (02) 4878 5592
PO Box 1167, Mittagong NSW 2575
enquiries@completemedicalaustralasia.com.au
www.completemedicalaustralasia.com.au
ABN: 63 169 647 787